

## MOSFET Patient Dose Verification Systems

Technology Built  
on 35 Years  
of Dosimetry  
Experience


• IMRT • IGRT/TOMOTHERAPY • VMAT • TBI • RADIOLOGY • BRACHYTHERAPY



### mobileMOSFET wireless FEATURES:

- Dose obtained in real-time
- Simplifies dosimetry and minimizes QA time
- Wireless with Bluetooth technology
- Export dosimetry data to R&V system

For more information  
and sales, please email to  
[bmcinfo@teambest.com](mailto:bmcinfo@teambest.com).



### Portable dosimeter FEATURES:

- Stand-alone dosimetry system
- Portable & compact with LCD display
- Optional data export via USB
- No software or computer necessary



### MOSFET dosimeter FEATURES:

- Isotropic (+2% for 360°)
- Active region of 0.2 x 0.2 mm
- Pinpoint measurement for high gradient fields
- Same calibration factor for photons and electrons at MeV energies